**SELECT DIGITAL BLOOD  
PRESSURE MONITOR**

Item # 04-630-001    Item # 04-631-001  
Item # 04-630-006    Item # 04-631-006

Please read this guidebook  
completely before operating the unit.

**DETAILED  
GUIDEBOOK**

*English • Spanish*

## LIMITED WARRANTY

---

The warrantor guarantees that its Digital Blood Pressure Monitor will be free from manufacturing defects under normal use for 5 years from the original purchase date of the product. This warranty covers only normal use, and does not apply for use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the Digital Blood Pressure Monitor. The warranty is voided if the Digital Blood Pressure Monitor product is misused or abused in any manner.

If the Digital Blood Pressure Monitor fails to operate during the first 5 years from the time originally purchased, please call our Customer Care Help Line at 800-622-4714 and we will try to resolve the issue as soon as possible. If the warrantor determines the unit failed to operate due to a manufacturing defect, the unit will be replaced at the option of the warrantor. The purchaser will be responsible for the prepaid postage along with \$9.99 to cover return shipping and insurance. The purchaser must return the original Blood Pressure Monitor prior to a replacement product being provided. Replacement of the product is the only remedy covered by the Limited 5-Year Warranty.

This warranty gives you specific legal rights, which vary from state to state. To be eligible for the warranty, please register your product online within 10 days of purchase date at [www.LiveHealthSmart.com](http://www.LiveHealthSmart.com).

We make no other express or implied warranties, arising by operation of law or otherwise, or any warranty of merchantability or fitness for a particular use or purpose whether or not the use or purpose has been disclosed to the warrantor in specifications, drawings or otherwise, and whether or not the warrantor's products are specifically designed and/or manufactured by the warrantor for the buyer's use or purposes, except for the limited warranty stated above. The warrantor will not be responsible for any indirect, incidental, special, consequential, or punitive damages or other loss, including, but not limited to, damage to or loss of other property or equipment and personal injuries, whether to purchaser or others. The warrantor shall in no event be liable to the purchaser for any amount in excess of the cost to replace the product.

# INDEX

---

Introduction . . . . .	3
Product Identification & Caution . . . . .	4
General Blood Pressure Information . . . . .	5-7
Important Information Before Use . . . . .	7
World Health Organization (WHO) Indicator . . . . .	8
Irregular Heartbeat Detection . . . . .	9
Battery Installation/Replacement . . . . .	10
Setting the Date/Time. . . . .	11
Applying Your Blood Pressure Cuff. . . . .	12
Taking Your Blood Pressure Reading. . . . .	13
Recalling Measurements in Memory . . . . .	14
Description of Symbols . . . . .	15
Troubleshooting . . . . .	16
Care and Maintenance . . . . .	17
Product Specifications . . . . .	18-25
Blood Pressure Logs . . . . .	26

## INTRODUCTION

---

Thank you for purchasing a Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.


To achieve the maximum benefit from your blood pressure monitor, we recommend that you first consult with your physician or trained healthcare professional.

The method of measurement that your Monitor uses is called the oscillometric method. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor easy to use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

## PRODUCT IDENTIFICATION

---

 **CAUTION:**

You can stop the inflation or deflation process anytime by pressing the POWER button.

# GENERAL BLOOD PRESSURE INFORMATION

---

## What is Blood Pressure?

Blood pressure is the pressure that is exerted by blood flowing against the walls of the blood vessels throughout your body.

Your heart, which is the center of the circulatory system, provides the force for the blood to flow or circulate. When your heart contracts or “beats,” the blood is forced through the blood vessels increasing the pressure. This is the highest pressure in the cycle or what is referred to as SYSTOLIC blood pressure. In-between beats your heart relaxes and your blood pressure decreases. This is referred to as DIASTOLIC blood pressure.

This complete series of events, which occurs in a single heartbeat, is known as the CARDIAC CYCLE.

Your Digital Blood Pressure Monitor will automatically read your blood pressure and display both systolic and diastolic readings on the screen. Your systolic will be positioned as the upper number and the diastolic reading will be the lower number.

(systolic) **120/80** (diastolic)

Blood pressure is measured in millimeters (mm) of mercury (Hg) and is generally recorded with the systolic pressure (120) listed first and the diastolic pressure (80) listed second. The numbers are typically separated by a slash mark (/) as shown.

Both pressure readings, the SYSTOLIC and DIASTOLIC, are necessary for a physician to evaluate the status of a patient’s blood pressure.

Please contact your physician for specific information regarding your own blood pressure.

## What Influences Blood Pressure?

Many factors such as genetics, age, sex, altitude, physical activity, anxiety, muscular development, certain medications, or even the time of day can influence blood pressure. Influences such as sleep or relaxation decrease blood pressure, while anxiety or exercise increase blood pressure.

## Why Monitor Your Blood Pressure at Home?

A visit to a physician's office can be stressful for a patient and anxiety is a known factor in raising blood pressure. This temporary occurrence of elevated blood pressure at the physician's office is commonly referred to as "white coat syndrome."

Whether or not you experience "white coat syndrome," home blood pressure monitoring provides you with the opportunity to supplement measurements taken at your physician's office. These home readings, when taken over a period of time, can show an accurate indication of change. Furthermore, your records can assist your physician in evaluating your health and in making important decisions in the diagnosis and treatment of your condition. Because of this, it is important to take consistent daily measurements of your blood pressure.

The variations in your individual readings should only be interpreted by your physician or a trained healthcare professional.

## WHO Blood Pressure Classifications


Standards for assessment of high or low blood pressure have been established by the World Health Organization (WHO) as shown on the following chart:

<b>Category</b>	<b>Systolic (mmHg)</b>	<b>Diastolic (mmHg)</b>	<b>Color Indicator</b>
Optimal	<120	<80	Green
Normal	120–129	80–84	Green
High - Normal	130–139	85-89	Green
Stage 1 - Mild			
Borderline Hypertension	140–159	90–99	Yellow
Stage 2 - Moderate			
Hypertension	160–179	100-109	Orange
Stage 3 - Severe Hypertension	≥180	≥110	Red
Isolated Systolic Hypertension	≥140	<90	
	140–159	<90	Yellow
	160–179	<90	Orange
	≥180	<90	Red

This chart is only a general guideline. Contact your physician or trained healthcare professional to determine your NORMAL blood pressure.

## Variations in Blood Pressure

Blood pressure is influenced by many factors and can change from moment to moment. Normally, blood pressure is lowest during sleeping periods and rises during the day. The graph below represents variations in blood pressure shown over a day with measurements taken every 5 minutes.

The dotted line represents the sleep period. The rise in blood pressure at 1 p.m. (A in the graph) reflects stressful occurrence, and at 7 p.m. (B in the graph), a period of exercise.


## IMPORTANT INFORMATION BEFORE USE

---

1. Blood pressure readings should be interpreted by a physician or trained healthcare professional who is familiar with your medical history.
2. Measure your blood pressure in a quiet place. You should be seated and relaxed.
3. Avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to measuring your blood pressure. If you are exhibiting signs of stress, avoid measuring your blood pressure until the feeling subsides.
4. Rest 15 minutes prior to measuring your blood pressure.
5. Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
6. Keep the monitor stable during measurements. Remain still; do not talk during measurements.
7. Record your daily readings on the chart provided.
8. Take your readings at the same time, each day (or as recommended by your physician).
9. Wait a minimum of 15 minutes between readings. The wait time may vary depending on individual physiological characteristics.
10. This device is intended for adult use.
11. The inflation or deflation process can be stopped by pressing the POWER button.

## WORLD HEALTH ORGANIZATION (WHO) INDICATOR

The WHO Indicator, located on the left side of the Display, will identify your reading within a range in accordance with the guidelines set by the World Health Organization. See the chart below for reference.

SYS  DIA 	<b>Stage 3 Hypertension</b> Systolic pressure: $\geq 180$ mmHg -or- Diastolic pressure: $\geq 110$ mmHg
SYS  DIA 	<b>Stage 2 Hypertension</b> Systolic pressure: 160-179 mmHg -or- Diastolic pressure: 100-109 mmHg
SYS  DIA 	<b>Stage 1 Hypertension</b> Systolic pressure: 140-159 mmHg -or- Diastolic pressure: 90-99 mmHg
SYS  DIA 	<b>High Normal</b> Systolic pressure: 130-139 mmHg -or- Diastolic pressure: 85-89 mmHg
SYS  DIA 	<b>Normal</b> Systolic pressure: 120-129 mmHg -or- Diastolic pressure: 80-84 mmHg
SYS  DIA 	<b>Optimal</b> Systolic pressure: $< 120$ mmHg -and- Diastolic pressure: $< 80$ mmHg


## IRREGULAR HEARTBEAT DETECTION

---

Your digital blood pressure monitor features Irregular Heartbeat Detection (IHB). This feature allows users to accurately monitor blood pressure even if an irregular heartbeat should occur.

When an irregular heartbeat is detected, the IHB icon will appear on the display, *Fig. 1*.

**NOTE:** Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat if this symbol appears frequently at the end of a reading.

*Fig. 1*


If Irregular Heartbeat caused by common arrhythmias is detected in the process of blood pressure measurement, a signal will be displayed. Under this condition, the Electronic Sphygmomanometers will function, but the results may not be accurate and it's suggested that you consult your physician for accurate assessment.

## BATTERY INSTALLATION/REPLACEMENT

1. Press down on tab and lift the cover.
2. Insert or replace 4 “AA” batteries into the compartment, matching the indicated polarity symbols, *Fig. 1*.
3. To close, match up and connect the bottom hooks of the cover, then push the top end.

The monitor, the batteries, and the cuff must be disposed of according to local regulations at the end of their usage.

*Fig. 1*


**NOTE:** Replace the batteries when the Low Battery symbol appears on the display or when the display remains blank when the POWER button is pressed.

4. It is recommended to remove the batteries if the unit will not be used for an extended period of time.

### AC Power Adapter (optional)

Use the AC adapter for battery-free operation. The AC adapter plugs into the side of the unit as shown, *Fig. 2*. Please remember to practice electrical safety precautions while using the adapter.

#### AC Adapter Port

## SETTING THE DATE/TIME

**NOTE:** When the unit is off, the date and time will still appear at the top of the screen.


1. When the unit is off, press and release both the POWER and MEMORY button at the same time. The month will begin to flash, *Fig. 1*.
2. Press the MEMORY button to change until the desired number appears.

**NOTE:** To scroll through the numbers more quickly, hold down the MEMORY button.


3. Press the POWER button to set the entry.
4. The date, *Fig. 2*, followed by the hour and minutes will flash. Repeat steps 2 and 3 to adjust the date/time setting.
5. After setting the minutes, the unit will automatically exit out of the date/time setting mode, *Fig. 3*.

**NOTE:** The time is expressed in 24-hour military time. The chart below explains the time conversion.


Regular Time (12 hours)	Military Time (24 hours)	Regular Time (12 hours)	Military Time (24 hours)
12 am Midnight	0	12 pm Noon	12
1 am	1	1 pm	13
2 am	2	2 pm	14
3 am	3	3 pm	15
4 am	4	4 pm	16
5 am	5	5 pm	17
6 am	6	6 pm	18
7 am	7	7 pm	19
8 am	8	8 pm	20
9 am	9	9 pm	21
10 am	10	10 pm	22
11 am	11	11 pm	23

*Fig. 1*

*Fig. 2*

*Fig. 3*

# APPLYING YOUR BLOOD PRESSURE CUFF

Before applying your blood pressure cuff, be sure you have selected the appropriate cuff:

**Standard Cuff** fits arm circumference 11¾" – 16½" (30 cm - 42 cm)

**Large Cuff** fits arm circumference 16½" – 18¾" (42 cm - 48 cm)

**NOTE:** If for any reason you are unable to or should not use your left arm, please modify the following instructions and apply the cuff to your right arm. Your physician can tell you which arm is best for you to use.


1. Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
2. Be seated with your feet flat on the floor.
3. Position the cuff on a solid surface with the tubing facing UP and AWAY from you, *Fig. 1*. The metal bar on the cuff should be to the left of the tubing.
4. Widen or open the cuff by pulling or rolling the bottom of the cuff towards the right, *Fig. 2*. This should open the cuff, creating a cylinder. Do not extend the cuff beyond the metal bar.
5. Insert your arm into the cuff (cylinder). Position the (Φ) mark over the main artery (on the inside of your arm), *Fig. 3*.
6. The bottom edge of the cuff should be positioned approximately one inch ABOVE the elbow joint, *Fig. 4*.
7. Reaching underneath your left arm with your right hand, pull the end of the cuff towards your body to tighten the cuff, *Fig. 5*. Wrap and secure the cuff making sure that the (Φ) mark remains as shown, *Fig. 3*.
8. The cuff should fit comfortably, yet snugly around your arm.
9. You should be able to insert two fingers between your arm and the cuff.


*Fig. 1*


*Fig. 2*

*Fig. 3*

*Fig. 4*


*Fig. 5*

**NOTE:** The cuff size is suitable for use when the vertical "Index" mark is within the horizontal "OK" range. A different sized cuff is needed if the "Index" mark is outside of the "OK" range.


# TAKING YOUR BLOOD PRESSURE READING

Proceed only after reading the previous sections of this manual.

1. Position the monitor on a flat, stable surface with the digital display panel in view. Insert the cuff tubing connector into the left side of the monitor, *Fig. 1*.
2. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, *Fig. 2*. Relax your left hand. Sit comfortably and don't cross your legs. The middle of the cuff should be at the level of the right atrium of the heart.
3. Press the POWER Button to turn the unit on.
4. This unit will run a self-test, *Fig. 3*. Select the memory bank (1 or 2) to store the reading in by pressing the MEMORY button. Press POWER to confirm your selection.

*Fig. 1*


*Fig. 2*


5. This unit will show the last measured reading, *Fig. 4*, or if there is no previous measurement, "0" as shown in *Fig. 5*.
6. When the display appears as shown, *Fig. 6*, the monitor is ready. The cuff will automatically begin to inflate.
7. Once the cuff inflates to the appropriate pressure level, the measurement will begin.
8. The unit will begin to automatically deflate and the numbers appearing on the display will decrease.
9. When the air is deflated, the measurement is complete. The unit will display your blood pressure (systolic and diastolic) measurements and pulse rate, *Fig. 7*. Record your readings in a chart.
10. Your measurements will automatically be saved in memory.
11. Press the POWER button to turn the power off, or the unit will automatically shut off after approximately 1 minute.
12. Disconnect the tubing from the monitor prior to storing.


*Fig. 3*

*Fig. 4*

*Fig. 5*

*Fig. 6*

*Fig. 7*

**NOTE:** This monitor inflates the cuff to approximately 190 mmHg or 40 mmHg higher than the last systolic pressure stored in memory. If the system detects that the cuff needs additional pressure, it will automatically re-inflate to approximately 40 mmHg higher than the last pressure level.


## RECALLING MEASUREMENTS IN MEMORY

You can recall up to 120 total measurements in memory (60 per memory bank). These measurements can be shared with your physician or trained healthcare professional.

1. Press and release the MEMORY button. The unit will display either U1 or U2 depending upon the user bank selected. It will also show how many measurements are stored within that bank.
2. To select the appropriate user bank press the POWER button to switch between U1 and U2.
3. Once the desired user memory bank is selected, press the MEMORY button again to scroll through the saved measurements, *Fig 1*. After you've scrolled through all the readings, the unit will automatically turn off.

*Fig. 1*



*Fig. 2*

**NOTE:** Each memory bank can store up to 60 readings. When the number of readings exceeds 60, the oldest data will be replaced with the new record.

### To erase all stored measurements:

1. Press and release the MEMORY button. The unit will display either U1 or U2 depending upon the user bank selected and how many measurements are stored within that bank.
2. Select the desired user bank by pressing and releasing the MEMORY button.
3. Press and hold down the MEMORY button until two horizontal lines flash, *Fig 2*. This indicates that all measurements within the selected memory bank have been erased.

## DESCRIPTION OF DISPLAY SYMBOLS

 <p>Indicates a measurement within Stage 3 Hypertension.</p>	 Shows the pulse rate per minute.
	 Occurs when an error was made during measurement.
 <p>Indicates a measurement within Stage 2 Hypertension.</p>	 Appears when there is an irregular heartbeat detected during measurement.
	 Appears when the cuff is unstable or it has too much air left in it. Once the arrow stops flashing, the unit is ready to take a measurement.
 <p>Indicates a measurement within Stage 1 Hypertension.</p>	 Appears when the unit is ready or a measurement error has occurred. Restart the measurement if an error occurs.
	 Appears when recalling a measurement from memory.
 Appears when batteries should be replaced.	

## DESCRIPTION OF SYMBOLS

Symbol for “THE OPERATION GUIDE MUST BE READ” (The sign background color: blue. The sign graphical symbol: white)

Symbol for “WARNING”

Symbol for “TYPE BF APPLIED PARTS” (The cuff is type BF applied part)

Symbol for “ENVIRONMENT PROTECTION” – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice.

SN

Symbol for “SERIAL NUMBER”

Symbol for “KEEP DRY”

## TROUBLESHOOTING

---

If any abnormality occurs during use, please check and correct the following:

<b>Condition</b>	<b>Correction</b>
Display is blank when power is on.	Check and correct the polarity of the installed batteries. Reinstall or replace batteries.
Measurement incomplete or abnormally low or high values displayed	Review and follow “Applying Your Blood Pressure Cuff” and “Taking Your Blood Pressure Reading” sections. Replace the batteries.
EE symbol displays.	Do not talk or move during the measurement. Review “Applying Your Blood Pressure Cuff” section.
Measurements are different from those typically measured by physician or every measurement reading is different.	Remember, blood pressure readings are influenced by physical and mental conditions and/or even the time of day. Daily measurements should be interpreted by your physician.

## CARE AND MAINTENANCE

---

1. Disconnect the cuff tubing from the monitor and remove batteries prior to storing.
2. Only use a soft, dry cloth to clean your blood pressure monitor.
3. Avoid using any types of liquids on the monitor or cuff.
4. Do not store the unit where it will be exposed to direct sunlight, dust or humidity. Avoid extreme temperatures.
5. Never disassemble the monitor or cuff.
6. Dropping or subjecting your blood pressure monitor to strong shocks should be avoided.
7. Remove batteries prior to storing.
8. It is recommended to disinfect the cuff twice a week if needed (for example, in hospital or in clinic). Wipe the inner side (the side that contacts skin) of the cuff with a soft cloth squeezed after moistened with Ethyl alcohol (75-90%), then dry the cuff by airing.
9. No component of the monitor should be serviced by the user. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information will assist the user's appropriately qualified technical personnel in repair. Those parts of equipment that are designated repairable can be supplied.
10. It is recommended for the performance to be checked every 2 years.
11. Please keep the cuff clean. If the cuff becomes dirty, remove it from the monitor and clean it by hand in a mild detergent, then rinse it thoroughly in cold water. Never dry the cuff in clothes dryer or iron it. Cleaning the cuff after the usage of every 200 times is recommended.
12. The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or 5 years.
13. The cuff's integrity is maintained after 1,000 open-close cycles of the closure.

# PRODUCT SPECIFICATIONS

---

<b>Name and item number</b>	Select Digital Blood Pressure Arm Monitor, 04-630-001, 04-630-006, 04-631-001, 04-631-006
<b>Display system</b>	Digital display/LCD
<b>Measuring method</b>	Oscillometric
<b>Power source</b>	--- 4 "AA" batteries or AC Adaptor
<b>Measuring range</b>	Pressure: 0-300 mmHg Pulse: 40-180 beats/minute
<b>Accuracy</b>	Pressure: $\pm 3$ mmHg Pulse: $\pm 5\%$ of reading
<b>Diastolic Range</b>	40 a 199 mmHG
<b>Systolic Range</b>	60 a 260 mmHG
<b>Inflation</b>	Automatic Inflation
<b>Deflation</b>	Automatic Pressure Release Valve/Manual fast release
<b>Memory</b>	Built-in memory enabling display of up to 120 measurements
<b>Automatic shut-off</b>	Approximately 1 minute after last button operation
<b>Operation environment</b>	Temperature 41°F – 104°F (5°C – 40°C) Humidity $\leq 90\%$
<b>Storage environment</b>	Temperature -4°F – 131°F (-20°C – 55°C) Humidity $\leq 95\%$
<b>Monitor dimensions</b>	5" x 5 $\frac{1}{8}$ " x 2 $\frac{3}{8}$ "
<b>Weight</b>	11.11 oz. (without batteries)
<b>Arm circumference range</b>	Adult cuff fits arm circumference: 11 $\frac{3}{4}$ " – 16 $\frac{1}{2}$ " (30 cm - 42 cm) Large cuff fits arm circumference 6 $\frac{1}{2}$ " – 18 $\frac{7}{8}$ " (42 cm - 48 cm)
<b>Accessories</b>	Cuff, batteries, case, detailed guidebook, quick ` guide, blood pressure log

Specifications are subject  
to change without notice.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

## ELECTROMAGNETIC COMPATIBILITY INFORMATION

**Table 1**

**For all ME EQUIPMENT and ME SYSTEMS**

Guidance and manufacture's declaration - electromagnetic emission		
The 04-630-001, 04-630-006, 04-631-001 and 04-631-006 are intended for use in the electromagnetic environment specified below. The customer or the user of the 04-630-001, 04-630-006, 04-631-001 and 04-631-006 should assure that it is used in such an environment		
Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The 04-630-001, 04-630-006, 04-631-001 and 04-631-006 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The 04-630-001, 04-630-006, 04-631-001 and 04-631-006 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3		

**Table 2**  
**For all ME EQUIPMENT and ME SYSTEMS**

Guidance and manufacture's declaration - electromagnetic immunity			
The 04-630-001, 04-630-006, 04-631-001 and 04-631-006 are intended for use in the electromagnetic environment specified below. The customer or the user of the 04-630-001 and 04-630-006 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance.
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>NOTE:</b> UT is the a.c. mains voltage prior to application of the test level.			

**Table 3**  
**For all ME EQUIPMENT and ME SYSTEMS (continued)**

Guidance and manufacture's declaration - electromagnetic immunity			
The 04-630-001, 04-630-006, 04-631-001 and 04-631-006 are intended for use in the electromagnetic environment specified below. The customer or the user of the 04-630-001, 04-630-006, 04-631-001 or 04-631-006 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test	Compliance level	Electromagnetic environment guidance.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the 04-630-001, 04-630-006, 04-631-001 or 04-631-006, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

**Table 3 - Continued**  
**For all ME EQUIPMENT and ME SYSTEMS**

<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V/m</p>	<p>Recommended separation distance:</p> <p><math>d=1.2 \sqrt{P}</math> 80 MHz to 800 MHz</p> <p><math>d=1.3 \sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level <sup>b</sup> in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: ((•))</p>
<p><b>NOTE 1:</b> At 80 MHz and 800 MHz, the higher frequency range applies.  <b>NOTE 2:</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

**Table 3 - Continued**  
**For all ME EQUIPMENT and ME SYSTEMS**

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 04-630-001, 04-630-006, 04-631-001, 04-631-006 are used exceeds the applicable RF compliance level above, the 04-630-001, 04-630-006, 04-631-001, 04-631-006 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the 04-630-001, 04-630-006, 04-631-001, 04-631-006.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Table 4**  
**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

Recommended separation distances between portable and mobile RF communications equipment and the 04-630-001, 04-630-006, 04-631-001 and 04-631-006.			
The 04-630-001, 04-630-006, 04-631-001 and 04-631-006 are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 04-630-001, 04-630-006, 04-631-001 or 04-631-006 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 04-630-001, 04-630-006, 04-631-001 or 04-631-006 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter.		
	150 kHz to 80 MHz $d=1,2 \sqrt{P}$	80 MHz to 800 MHz $d=1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d=1,2 \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

**Table 4 - Continued**

**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people and animals.

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