Health♥Smart R_x

*Automatic Digital
Blood Pressure Monitor*

04-525-001

Please read this instruction manual
completely before operating this unit.

Limited One-Year Warranty

The warrantor guarantees that its Digital Blood Pressure Monitor will be free from manufacturing defects under normal use for a period of one year from the date of purchase subject to the following terms and conditions.

This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the Digital Blood Pressure Monitor. This warranty is voided if the Digital Blood Pressure Monitor product is misused or abused in any manner.

If the Digital Blood Pressure Monitor fails to operate during the time the original purchaser owns it, return it postage prepaid along with \$5.00 for return shipping and insurance to: Repair Department, 1931 Norman Drive South, Waukegan, IL 60085. When returning a product, please also include a letter with your name, address, phone number and brief description of the specific problem. Please carefully package the product to avoid any damage that may occur while in the transit. Packaging insurance with returned receipt is recommended. The warrantor will repair or replace the defective unit at the option of the warrantor. Repair or replacement of the defective unit, at the warrantor's option, is the sole remedy under this warranty.

ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO THE TIME THAT THE ORIGINAL CONSUMER PURCHASER OWNS THE PRODUCT. Some states do not allow limitations or how long an implied warranty lasts, so the above limitation may not apply to you.

This warranty constitutes the warrantor's only responsibility and obligation to repair and/or replace materials or components, or refund the purchase price. **The warrantor will not be responsible for any indirect, incidental, special, consequential, or punitive damages or other loss**, including, but not limited to, damage to or loss of other property or equipment and personal injuries, whether to purchaser or others. The warrantor shall in no event be liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchase price of the unit. Some states do not allow the exclusive or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

Toll-Free Customer Care Help Line: 1-800-622-4714
Monday - Friday 8:00 am - 4:30 pm CST

HealthSmart RX
1931 Norman Drive South • Waukegan, IL 60085
847-680-6811 • FAX 847-680-9646

Index

Introduction	3
Product Identification & Caution	4
Why Monitor Your Blood Pressure at Home?	5
General Blood Pressure Information	5-7
Important Information Before Use	7
Battery Installation/Replacement	8
Applying Your Blood Pressure Cuff	9
Taking Your Blood Pressure Reading	10
Description of Display Symbols	11
Troubleshooting	11
Care and Maintenance	12
Product Specifications	12
Sample Blood Pressure Log	13-14

Introduction


To achieve the maximum benefit from your blood pressure monitor, we recommend that you first consult with your physician or trained healthcare professional.

Thank you for purchasing a HealthSmart Automatic Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

The method of measurement that your Automatic Monitor uses is called the oscillometric method. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor easy to use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Product Identification and Caution

CAUTION:

You can stop the inflation or deflation process anytime by pressing the On/Off button.

Why Monitor Your Blood Pressure at Home

Many people experience "White Coat Syndrome", which is a temporary elevation of blood pressure during a visit at their physician's office. Although only temporary, it can cause blood pressure to rise and give a false indication of a person's true blood pressure.

Whether or not you experience "White Coat Syndrome", home blood pressure monitoring provides:

- 1) The opportunity to supplement your physician's office measurements.
- 2) An accurate indication of any change.
- 3) A record to assist your physician in evaluating your health and making important decisions.
- 4) Consistent, daily measurements of your blood pressure.

NOTE: Variations in your individual readings should only be interpreted by your physician or trained healthcare professional.

General Blood Pressure Information

What is Blood Pressure?

Blood pressure is the pressure that is exerted by blood flowing against the walls of the blood vessels throughout your body.

Your heart, which is the center of the circulatory system, provides the force for the blood to flow or circulate. When your heart contracts or beats, the blood is forced through the blood vessels increasing the pressure. This is the highest pressure in the cycle or what is referred to as SYSTOLIC blood pressure. In between beats, your heart relaxes and your blood pressure decreases. This is referred to as DIASTOLIC blood pressure.

This complete series of events, which occurs in a single heartbeat, is known as the CARDIAC CYCLE.

Your MABIS Automatic Blood Pressure Monitor will automatically read your blood pressure and display both systolic and diastolic readings on the screen. Your systolic will be positioned as the upper number and the diastolic reading will be the lower number.

(systolic)
120/80
(diastolic)

Blood pressure is measured in millimeters (mm) of mercury (Hg) and is generally recorded with the systolic pressure (120) listed first and the diastolic pressure (80) listed second. The numbers are typically separated by a slash mark (/) as shown.

Both pressure readings, the SYSTOLIC and DIASTOLIC, are necessary for a physician to evaluate the status of a patient's blood pressure.

Please contact your physician for specific information regarding your own blood pressure.

What Influences Blood Pressure?

Many factors such as genetics, age, sex, altitude, physical activity, anxiety, muscular development, certain medications or even the time of day can influence blood pressure. Influences such as sleep or relaxation decrease blood pressure, while anxiety or exercise increase blood pressure.

WHO Blood Pressure Classifications


Standards for assessment of high or low blood pressure have been established by the World Health Organization (WHO) as shown on the following chart:

CATEGORY	SYSTOLIC (mmHg)	DIASTOLIC (mmHg)
Optimal	<120	<80
Normal	120-129	80-84
High - Normal	130-139	85-89
Mild Hypertension	140-159	90-99
Moderate Hypertension	160-179	100-109
Severe Hypertension	≥180	≥110
Isolated Systolic Hypertension	≥140	<90

This chart is only a general guideline. Contact your physician or trained healthcare professional to determine your NORMAL blood pressure.

Variations in Blood Pressure

Blood pressure is influenced by many factors and can change from moment to moment. Normally, blood pressure is lowest during sleeping periods and rises during the day. The graph below represents variations in blood pressure shown over a day with measurements taken every 5 minutes.

The dotted line represents the sleep period. The rise in blood pressure at 1 p.m. (A in the graph) corresponds to a stressful occurrence and at 7 p.m. (B in the graph) a period of exercise.

Important Information Before Use

- Blood pressure readings should be interpreted by a physician or trained healthcare professional, who is familiar with your medical history.
- Perform your measurement in a quiet place. You should be seated and relaxed.
- Avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading. If you are exhibiting signs of stress, avoid taking your measurement until the feeling subsides.
- Rest 15 minutes prior to taking a reading.
- Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
- Keep the monitor stable during measurements. Remain still; do not talk during measurements.
- Record your daily readings on the chart provided.
- Take your readings at the same time, each day (or as recommended by your physician).
- Wait a minimum of 15 minutes between readings. The wait time may vary depending on individual physiological characteristics.
- This device is intended for adult use.
- The inflation or deflation process can be stopped by pressing the On/Off button.

Battery Installation/Replacement

- 1) Slide battery cover off, *Fig. 1*.
- 2) Insert or replace 4 “AA” batteries into the compartment, matching the indicated polarity symbols, *Fig. 1*.
- 3) Replace the battery cover.

NOTE: Replace the batteries when the Low Battery symbol appears on the display or when the display remains blank when the On/Off button is pressed.

- 4) It is recommended to remove the batteries if the unit will not be used for an extended period of time.


Fig. 1

AC Power Adapter (optional)

An AC adapter is available separately for battery-free operation. The AC adapter plugs into the side of the unit as shown, *Fig. 2*. Please remember to practice electrical safety precautions while using the adapter. To order an AC adapter, please see the enclosed order form.

AC Adapter

Output: 6V DC 600mA

AC Adapter Port

Fig. 2

Applying Your Blood Pressure Cuff

Before applying your blood pressure cuff, be sure you have selected the appropriate cuff:

Small Adult Cuff fits arm circumference 8-5/8” - 11-3/4” (22-30 cm)

Adult Cuff fits arm circumference 11-3/4” - 16-1/2” (30-42 cm)

Large Adult Cuff fits arm circumference 16-1/2” - 18-7/8” (42-48 cm)

NOTE: If for any reason you are unable to or should not use your left arm, please modify the following instructions and apply the cuff to your right arm. Your physician can tell you which arm is best for you to use.

1. Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
2. Be seated with your feet flat on the floor.
3. Position the cuff on a solid surface with the tubing facing UP and AWAY from you, *Fig. 1*. The metal bar on the cuff should be to the left of the tubing.

Fig. 1

4. Widen or open the cuff by pulling or rolling the bottom of the cuff towards the right, *Fig. 2*. This should open the cuff, creating a cylinder. Do not extend the cuff beyond the metal bar.
5. Insert your arm into the cuff (cylinder). Position the (Φ) mark over the main artery (on the inside of your arm), *Fig. 3*.
6. The bottom edge of the cuff should be positioned approximately one inch ABOVE the elbow joint, *Fig. 4*.

Fig. 2


Fig. 3

7. Reaching underneath your left arm with your right hand, pull the end of the cuff towards your body to tighten the cuff, *Fig. 5*. Wrap and secure the cuff making sure that the (Φ) mark remains as shown, *Fig. 3*.


Fig. 4

8. The cuff should fit comfortably, yet snugly around your arm. You should be able to insert two fingers between your arm and the cuff.




Fig. 5

NOTE: The cuff size is suitable for use when the vertical “Index” mark is within the horizontal “OK” range. A different sized cuff is needed if the “Index” mark is outside of the “OK” range.

Taking Your Blood Pressure Reading

Proceed only after reading the previous sections of this manual.

- Position the monitor on a flat, stable surface with the digital display panel in view. Insert the cuff tubing connector into the left side of the monitor, *Fig. 1*.
- Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, *Fig. 2*. Relax your left hand.
- Press the On/Off button to turn the unit on.
- This unit will run a self-test, *Fig. 3*.
- When the display appears as shown, *Fig. 4*, the monitor is ready.
- Press the Start button and the monitor will automatically begin to inflate.
- The display will show the increasing pressure in the cuff.
- The cuff will automatically inflate to approximately 187 mmHg.
- Once the cuff inflates to approximately 187 mmHg, the measurement will begin.
- The unit will begin to automatically deflate and the numbers appearing on the display will decrease.
- When all the air is deflated, the measurement is complete. The unit will display your blood pressure (systolic and diastolic) measurements and pulse rate, *Fig. 5*.
- Record you readings on a chart.
- To conserve the batteries, press the On/Off button to turn the power off. Otherwise, the unit will automatically shut off after approximately 3 minutes.
- Disconnect the tubing from the monitor prior to storing.

NOTE: This monitor inflates the cuff to approximately 187 mmHg. If the system detects that the cuff needs additional pressure, it will automatically reinflate to approximately 225 mmHg, 263 mmHg, then finally if needed 300 mmHg.

Description of Display Symbols

	Appears when batteries should be replaced.
	Shows the pulse rate per minute.
	Occurs when an error was made during measurement.
	Appears when the unit is ready or a measurement error has occurred. Restart the measurement if an error occurs.
	Appears when the cuff is unstable or it has too much air left in it. Once the arrow stops flashing, the unit is ready to take a measurement.
	Appears when the measurement is being taken.

Troubleshooting

If any abnormality occurs during use, please check and correct the following:

Condition	Correction
Display is blank when power is on.	Check and correct the polarity of the installed batteries. Reinstall or replace batteries.
Measurement incomplete or abnormally low or high values displayed.	Review and follow "Applying Your and "Taking Your Blood Pressure Reading" sections. Replace the batteries.
EE and RETEST symbol displays.	More pressure is needed in the cuff. Retake your measurement.
	Do not talk or move during the measurement. Review "Applying Your Blood Pressure Cuff" section.
	Pressure quickly fell. Check to make sure the cuff is correctly connected and that no air is leaking.
Measurements are different from those typically measured by physician or every measurement reading is different.	Remember, blood pressure readings are influenced by physical and mental conditions and/or even the time of day. Daily measurements should be interpreted by your physician.

Care and Maintenance

1. Disconnect the cuff tubing from the monitor and remove batteries prior to storing.
2. Only use a soft, dry cloth to clean your blood pressure monitor.
3. Avoid using any types of liquids on the monitor or cuff.
4. Do not store the unit where it will be exposed to direct sunlight, dust or humidity.
5. Avoid extreme temperatures.
6. Never disassemble the monitor or cuff.
7. Dropping or subjecting your blood pressure monitor to strong shocks should be avoided.

Product Specifications

Name and model number	Automatic Digital Blood Pressure Arm Monitor, 04-525-001
Display system	Digital display/LCD
Measuring method	Oscillometric
Power source	4 "AA" batteries
Measuring range	Pressure: 0-300 mmHg Pulse: 30-180 beats/minute
Accuracy	Pressure: ±3 mmHg Pulse: ±5% of reading
Inflation	Pump Driven
Deflation	Automatic Pressure Release Valve
Automatic shut-off	Approximately 3 minutes after last button operation
Operation environment	Temperature 41°F – 104°F (5°C – 40°C) Humidity < 85%
Storage environment	Temperature -4°F – 149°F (-20°C – 65°C) Humidity < 90%
Monitor dimensions	5-5/8" x 4" x 2-1/2"
Weight	13.4 oz. (without batteries)
Arm circumference range	Adult size cuff fits arm circumference: 11-3/4" - 16-1/2" (30-42 cm),
Accessories	Cuff, detailed guidebook, quick start guide, blood pressure log, four "AA" batteries, storage case
Options (sold separately)	Large adult size cuff fits arm circumference: 16-1/4" -18-7/8" (42-48 cm), Small adult size cuff fits arm circumference: 8-5/8" -11-3/4" (22-30 cm), AC Adapter

Specifications are subject to change without notice.

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©2008 HealthSmart RX