healthsmart[®]

**PREMIUM TALKING DIGITAL WRIST
BLOOD PRESSURE MONITOR**

Item # 04-895-001

Please read this guidebook
completely before operating the unit.

**DETAILED
GUIDEBOOK**

English • Spanish

LIMITED 5-YEAR WARRANTY

The warrantor guarantees that its Digital Blood Pressure Monitor will be free from manufacturing defects under normal use for 5 years from the original purchase date of the product. This warranty covers only normal use, and does not apply for use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the Digital Blood Pressure Monitor. The warranty is voided if the Digital Blood Pressure Monitor product is misused or abused in any manner.

If the Digital Blood Pressure Monitor fails to operate during the first 5 years from the time originally purchased, please call our Customer Care Help Line at 800-622-4714 and we will try to resolve the issue as soon as possible. If the warrantor determines the unit failed to operate due to a manufacturing defect, the unit will be replaced at the option of the warrantor. The purchaser will be responsible for the prepaid postage along with \$9.99 to cover return shipping and insurance. The purchaser must return the original Blood Pressure Monitor prior to a replacement product being provided. Replacement of the product is the only remedy covered by the Limited 5-Year Warranty.

This warranty gives you specific legal rights, which vary from state to state. To be eligible for the warranty, please register your product online within 10 days of purchase date at www.LiveHealthSmart.com.

We make no other express or implied warranties, arising by operation of law or otherwise, or any warranty of merchantability or fitness for a particular use or purpose whether or not the use or purpose has been disclosed to the warrantor in specifications, drawings or otherwise, and whether or not the warrantor's products are specifically designed and/or manufactured by the warrantor for the buyer's use or purposes, except for the limited warranty stated above. The warrantor will not be responsible for any indirect, incidental, special, consequential, or punitive damages or other loss, including, but not limited to, damage to or loss of other property or equipment and personal injuries, whether to purchaser or others. The warrantor shall in no event be liable to the purchaser for any amount in excess of the cost to replace the product.

INDEX

Introduction	2
Product Identification & Caution	3
General Blood Pressure Information	4-6
Important Information Before Use	7
World Health Organization (WHO) Indicator	8
Irregular Heartbeat Detection	9
Battery Installation/Replacement	10
Setting the Date/Time	11
Voice Language Setting.	12-13
Applying Your Blood Pressure Cuff	14-15
Taking Your Blood Pressure Reading.	16-17
Recalling Measurements in Memory	18
Description of Display Symbols	19-21
Troubleshooting	22
Care and Maintenance	23
Product Specifications	24-31
Sample Blood Pressure Logs	32
Blood Pressure Log	33

INTRODUCTION

Thank you for purchasing a Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

To achieve the maximum benefit from your blood pressure monitor, we recommend that you first consult with your physician or trained healthcare professional.


Your Digital Monitor uses oscillometric method of measurement. The monitor detects your blood's movement through the artery in your wrist and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor easy-to-use. The cuff is preformed for the wrists and is easy to apply.

Clinical research has proven a direct relationship between blood pressure in the wrist and in the arm. Because the arteries in the wrist and arm are connected, changes in wrist blood pressure reflect changes in arm blood pressure. Frequent wrist blood pressure measurements can provide you and your physician with an accurate indication of change in your true blood pressure.

People with severe hypertension, severe arteriosclerosis or severe diabetes should consult their physician about monitoring blood pressure at the wrist.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

PRODUCT IDENTIFICATION AND CAUTION

CAUTION

While taking a measurement, you can stop the inflation or deflation process of the cuff at any time by pressing the POWER button.

GENERAL BLOOD PRESSURE INFORMATION

What is Blood Pressure?

Blood pressure is the pressure that is exerted by blood flowing against the walls of the blood vessels throughout your body.

Your heart, which is the center of the circulatory system, provides the force for the blood to flow or circulate. When your heart contracts or beats, the blood is forced through the blood vessels increasing the pressure. This is the highest pressure in the cycle or what is referred to as **SYSTOLIC BLOOD PRESSURE**. In between beats your heart relaxes and your blood pressure decreases. This is referred to as **DIASTOLIC BLOOD PRESSURE**.

This complete series of events which occur in a single heartbeat is known as the **CARDIAC CYCLE**.

Your Digital Blood Pressure Monitor will automatically read your blood pressure and display both systolic and diastolic readings on the screen. Your systolic will be positioned as the upper number and the diastolic reading will be the lower number.

(systolic) **120/80** (diastolic)

Blood pressure is measured in millimeters (mm) of Mercury (Hg) and is generally recorded with the systolic pressure (120) listed first and the diastolic pressure (80) listed second. The numbers are typically separated by a slash mark (/) as shown above.

Both pressure readings, the **SYSTOLIC** and the **DIASTOLIC**, are necessary for a physician to evaluate the status of a patient's blood pressure.

Please contact your physician for specific information regarding your own blood pressure.

What Influences Blood Pressure?

Many factors such as genetics, age, sex, altitude, physical activity, anxiety, muscular development, certain medications, or even the time of day can influence blood pressure. Influences such as sleep or relaxation decrease blood pressure, while anxiety or exercise increases blood pressure.

Why Monitor Your Blood Pressure at Home?

A visit to a physician's office can be stressful for the patient and anxiety is a known factor in raising blood pressure. This temporary occurrence of elevated blood pressure at the physician's office is commonly referred to as "white coat syndrome."

Whether or not you experience "white coat syndrome," home blood pressure monitoring provides you with the opportunity to supplement measurements taken at your physician's office. These home readings, when taken over a period of time, can show an accurate indication of change. Furthermore, your records can assist your physician in evaluating your health and in making important decisions in the diagnosis and treatment of your condition. Because of this, it is important to take consistent daily measurements of your blood pressure.

The variations in your individual readings should only be interpreted by your physician or a trained healthcare professional.

WHO Blood Pressure Classifications

Standards for assessment of high or low blood pressure, regardless of age, have been established by the World Health Organization (WHO) as shown in the chart on the following page.


GENERAL BLOOD PRESSURE INFORMATION (continued)

Category	Systolic (mmHg)	Diastolic (mmHg)	Color Indicator
Optimal	<120	<80	Green
Normal	120–129	80–84	Green
High - Normal	130–139	85–89	Green
Stage 1 - Mild			
Borderline Hypertension	140–159	90–99	Yellow
Stage 2 - Moderate			
Hypertension	160–179	100–109	Orange
Stage 3 - Severe Hypertension	≥180	≥110	Red
Isolated Systolic Hypertension	≥140	<90	
	140–159	<90	Yellow
	160–179	<90	Orange
	≥180	<90	Red

This chart is only a general guideline. Contact your physician to determine your NORMAL blood pressure.

Variations in Blood Pressure

Blood pressure can be influenced by many factors. This change can occur from moment to moment. Normally, blood pressure is lowest during the sleeping period and rises during the day. The graph below represents the variations in blood pressure shown over a day with measurements taken every 5 minutes.

The dotted line represents the sleep period. The rise in blood pressure at 1 p.m. (A in the graph) reflects stressful occurrence, and at 7 p.m. (B in the graph), a period of exercise.

IMPORTANT INFORMATION BEFORE USE


1. Blood pressure measurements should only be interpreted by a physician or a trained healthcare professional who is familiar with your medical history. Through regular use of this device and recording of your measurements, you can keep your physician informed of the changes in your blood pressure.
2. Perform your measurement in a quiet place. You should be seated in a relaxed position.
3. Avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading. If you are exhibiting signs of stress, avoid taking your measurement until the feeling subsides.
4. Rest 15 minutes prior to taking a reading.
5. Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
6. Apply the cuff to the proper position on the left wrist. The unit should be located on the inside of your left wrist with the display panel in view.
7. Keep the monitor stable during measurement to achieve an accurate reading. Remain still; do not talk during the measurement.
8. Record your daily readings on the chart provided.
9. Take your readings at the same time each day (or as recommended by your physician) to get an accurate indication of change in your true blood pressure.
10. Wait a minimum of 15 minutes between readings to allow for the blood vessels to return to normal. The wait time may vary depending on your individual physiological characteristics.
11. Although such cases are rare, for those with an extremely weak pulse or irregular pulse, errors may result which prevent proper measurement. If abnormal variations are noticed, consult with your physician or trained healthcare professional.
12. This device is intended for adult use.
13. It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.
14. Please do not use the cuff other than supplied by the manufacturer, otherwise it may bring biocompatible hazard and might result in measurement error.
15. Please do not share the cuff with other infective person to avoid cross-infection.

16. This Electronic Sphygmomanometers should never be used on infants or young children. Consult your physician or other health care professionals before use on older children.
17. Prolonged over-inflation (cuff pressure exceed 300 mmHg or maintained above 15 mmHg for longer than 3 minutes) of the bladder may cause ecchymoma of your arm.
18. The application of the cuff over a wound or inflammation diseases.
19. The application of the cuff on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present;
20. The application of the cuff on the arm on the side of a mastectomy;
21. Simultaneously used with other monitoring medical equipments on the same limb.
22. Need to check the blood circulation of the user.
23. The monitor might not meet its performance specifications or cause safety hazard if stored or used outside the specified temperature and humidity ranges in specifications.
24. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.

While taking a measurement, you can stop the inflation or deflation process of the cuff at any time by pressing the POWER button.

WORLD HEALTH ORGANIZATION (WHO) INDICATOR

The WHO Indicator, located on the left side of the Display, will identify your reading within a range in accordance with the guidelines set by the World Health Organization. See the chart below for reference.

SYS DIA	 Stage 3 Hypertension Systolic pressure: ≥ 180 mmHg -or- Diastolic pressure: ≥ 110 mmHg
SYS DIA	 Stage 2 Hypertension Systolic pressure: 160-179 mmHg -or- Diastolic pressure: 100-109 mmHg
SYS DIA	 Stage 1 Hypertension Systolic pressure: 140-159 mmHg -or- Diastolic pressure: 90-99 mmHg
SYS DIA	 High Normal Systolic pressure: 130-139 mmHg -or- Diastolic pressure: 85-89 mmHg
SYS DIA	 Normal Systolic pressure: 120-129 mmHg -or- Diastolic pressure: 80-84 mmHg
SYS DIA	 Optimal Systolic pressure: < 120 mmHg -and- Diastolic pressure: < 80 mmHg

IRREGULAR HEARTBEAT DETECTION

Your digital blood pressure monitor features Irregular Heartbeat Detection (IHB). This feature allows users to accurately monitor blood pressure even if an irregular heartbeat should occur.

When an irregular heartbeat is detected, the IHB icon will appear on the display, *Fig. 1*.

NOTE: Please consult your physician or trained healthcare professional for further information regarding an irregular heartbeat if this symbol appears frequently at the end of a reading.


Fig. 1

If Irregular Heartbeat caused by common arrhythmias is detected in the process of blood pressure measurement, a signal will be displayed. Under this condition, the Electronic Sphygmomanometers will function, but the results may not be accurate and it's suggested that you consult your physician for accurate assessment.


BATTERY INSTALLATION/REPLACEMENT

Replace the batteries when the Low Battery symbol appears on the display, *Fig. 1*, or when the display remains blank when the POWER Button is pressed.

To change the batteries:

1. Open the battery cover, *Fig. 2*.
2. Making sure the battery polarity is correct, install two AAA alkaline batteries, *Fig. 3*.
3. Close the battery cover, *Fig. 4*.

Fig. 1

The monitor, the batteries and the cuff must be disposed of according to local regulations at the end of their usage.

NOTE: It is recommended to remove the batteries if the unit will not be used for an extended period of time.


Fig. 2


Fig. 3


Fig. 4

SETTING THE DATE/TIME

NOTE: When the unit is off, the date and time will still appear at the top of the screen.

1. When the unit is off, press and release both the POWER and MEMORY button at the same time. The month will begin to flash, *Fig. 1*.
2. Press the MEMORY button to change until the desired number appears.

NOTE: To scroll through the numbers more quickly, hold down the MEMORY button.

3. Press the POWER button to set the entry.
4. The date, *Fig. 2*, followed by the hour and minutes will flash. Repeat steps 2 and 3 to adjust the date/time setting.
5. After setting the minutes, the unit will automatically exit out of the date/time setting mode, *Fig. 3*.

NOTE: The time is expressed in 24-hour military time. The chart below explains the time conversion.

Regular Time (12 hours)	Military Time (24 hours)	Regular Time (12 hours)	Military Time (24 hours)
12 am Midnight	0	12 pm Noon	12
1 am	1	1 pm	13
2 am	2	2 pm	14
3 am	3	3 pm	15
4 am	4	4 pm	16
5 am	5	5 pm	17
6 am	6	6 pm	18
7 am	7	7 pm	19
8 am	8	8 pm	20
9 am	9	9 pm	21
10 am	10	10 pm	22
11 am	11	11 pm	23


Fig. 1


Fig. 2



Fig. 3


VOICE LANGUAGE SETTING

To set the desired voice language:

1. When the unit is off, you can select the desired voice language by pressing and holding the MEMORY button.
2. Now the screen will cycle through L0, L1 and L2 as shown in *Fig. 1*.
L0 - voice function is off
L1 - voice is in English
L2 - voice is in Spanish
The monitor will say the corresponding language with the displayed language code.
3. Once the desired language has been selected, release the MEMORY button.

4. After the desired language has been selected, you can adjust the volume by pressing and holding down the MEMORY button until the bars on the side of the screen begin to increase *Fig. 2*.
5. Continue to hold the MEMORY button down until the volume reaches the desired level, then release the MEMORY button to set.
6. The unit will automatically turn off after the desired volume level has been set.

APPLYING YOUR BLOOD PRESSURE CUFF

Cuff Dimensions	3¼" x 2¾" x 1⅞" (8.26 cm x 7 cm x 2.9 cm)
Wrist Circumference Range	Suitable for 5½" – 7⅝" (14 cm – 19.5 cm)

Avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to measuring your blood pressure. If you are exhibiting signs of stress, avoid taking your measurement until the feeling subsides. Rest 15 minutes prior to taking a reading.

1. Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
2. You should be seated with your feet flat on the floor.
3. Hold your left arm in front of you with your palm facing upward.
4. Apply the preformed cuff to your left wrist. The digital display panel should be positioned on the inside area of your wrist facing you, *Fig. 1*.
5. Adjust the cuff as shown in *Fig. 2* from the edge of the 'Head of the Ulna Bone.'


Fig. 1


Fig. 2

6. Wrap and fasten the cuff securely yet comfortably around your wrist, *Fig. 3*.

NOTE: The left wrist should be used unless this is not a viable option for you. In consultation with your physician, you may modify the following instructions and apply the cuff to your right wrist.




Fig. 3

TAKING YOUR BLOOD PRESSURE READING

Proceed only after reading the previous sections of this manual.

1. Rest your elbow on a solid surface with your palm facing upward. **Elevate your hand so that the cuff is at the same level as your heart**, Fig. 1. Relax your left hand. Sit comfortably and don't cross your legs. The middle of the cuff should be at the level of the right atrium of the heart.
2. Press the POWER button to turn the unit on and select the memory bank you want your reading to be stored in by pressing the MEMORY button to alternate between memory bank 1 or 2. Press the POWER button to confirm your selection.
3. The monitor will begin to guide you with voice prompts if the voice function is turned on. This unit will run a self-test, Fig. 2, and show the last measured reading, Fig. 3. Or if there is no previous measurement, a "0" will be displayed as shown in Fig. 4.
4. When the display appears as shown, Fig. 5, the monitor is ready and will begin to guide you with voice prompts if the voice function is on. The cuff will automatically begin to inflate.
5. The display will show the increasing pressure in the cuff along with the measurement progress on the left side of the screen.
6. The cuff will automatically inflate to approximately 190 mmHg or 40 mmHg higher than the last systolic pressure stored in memory.
7. Once the cuff inflates to the appropriate pressure level, the measurement will begin.

- The unit will begin to automatically deflate and the numbers appearing on the display will decrease.
- When all the air is deflated, the measurement is complete. The unit will display your blood pressure (systolic and diastolic) measurements and pulse rate, *Fig. 6*.

With the voice function on (L1 or L2), the monitor will say the measurement with the language you selected according to the classification standard of WHO (World Health Organization).


Fig. 6

- Record your readings on a chart.
- At this point, your measurements will be automatically saved in memory.
- To conserve the batteries, press the POWER button to turn the power off. Otherwise, the unit will automatically shut off after approximately 1 minute.

NOTE: This monitor inflates the cuff to approximately 190 mmHg or 40 mmHg higher than the last systolic pressure stored in memory. If the system detects that the cuff needs additional pressure, it will automatically reinflate to approximately 40 mmHg higher than the last pressure level.

RECALLING MEASUREMENTS IN MEMORY

You can recall up to 120 total measurements in memory (60 per memory bank). These measurements can be shared with your physician or trained healthcare professional.

1. Press and release the MEMORY button. The unit will display either U1 or U2 depending upon the user bank selected. It will also show how many measurements are stored within that bank, *Fig. 1*.
2. To select the appropriate user bank press the POWER button to switch between U1 and U2.
3. Once the desired user memory bank is selected, press the MEMORY button again to scroll through the saved measurements. After you've scrolled through all the readings, the unit will automatically turn off.

NOTE: Each memory bank can store up to 60 readings. When the number of readings exceeds 60, the oldest data will be replaced with the new record.

To erase all stored measurements:

1. Press and release the MEMORY button. The unit will display either U1 or U2 depending upon the user bank selected and how many measurements are stored within that bank.
2. Select the desired user bank by pressing and releasing the MEMORY button.
3. Press and hold down the MEMORY button until two horizontal lines flash. This indicates that all measurements within the selected memory bank have been erased.


Fig. 1




Fig. 2


DESCRIPTION OF DISPLAY SYMBOLS

Indicates a measurement within Stage 3 Hypertension.

Indicates a measurement within Stage 2 Hypertension.

Indicates a measurement within Stage 1 Hypertension.

Appears when batteries should be replaced.

DESCRIPTION OF DISPLAY SYMBOLS (continued)

Shows the pulse rate per minute.

Occurs when an error was made during measurement.

Appears when there is an irregular heartbeat detected during measurement.

Appears when the cuff is unstable or it has too much air left in it. Once the arrow stops flashing, the unit is ready to take a measurement.

Appears when a measurement error has occurred. Restart the measurement.

Appears when recalling a measurement from memory.

Indicates the memory user's memory bank.

DESCRIPTION OF SYMBOLS

Symbol for “THE OPERATION GUIDE MUST BE READ” (The sign background color: blue. The sign graphical symbol: white)

Symbol for “WARNING”

Symbol for “TYPE BF APPLIED PARTS” (The cuff is type BF applied part)

Symbol for “ENVIRONMENT PROTECTION” – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice.

SN

Symbol for “SERIAL NUMBER”

Symbol for “KEEP DRY”

TROUBLESHOOTING

If any abnormality occurs during use, please check and correct the following:

Condition	Correction
Display is blank when power is on.	Check and correct the polarity of the installed batteries. Reinstall or replace batteries.
Measurement incomplete or abnormally low or high values displayed.	Review and follow “Applying your Blood Pressure Cuff” and “Taking Your Blood Pressure Reading” sections. Replace the batteries.
Measurements are different from those typically measured by physician or every measurement reading is different.	Remember, blood pressure readings are influenced by physical and mental conditions and/or even the time of day. Daily measurement recordings should be interpreted by your physician.
EE Symbol Displays.	Do not talk or move during the measurement. Review “Applying Your Blood Pressure Cuff” section.

CARE AND MAINTENANCE

1. Only use a soft, dry cloth to clean your blood pressure monitor.
2. Avoid using any types of liquids on the monitor or cuff.
3. Store your monitor in the provided storage case when not in use.
4. Do not store the unit where it will be exposed to direct sunlight, dust or humidity. Avoid extreme temperatures.
5. Never disassemble the monitor or cuff.
6. Dropping or subjecting your blood pressure monitor to strong shocks should be avoided.
7. Remove batteries prior to storing.
8. It is recommended to disinfect the cuff twice a week if needed (for example, in hospital or in clinic). Wipe the inner side (the side that contacts skin) of the cuff with a soft cloth squeezed after moistened with Ethyl alcohol (75-90%), then dry the cuff by airing.
9. No component of the monitor should be serviced by the user. The circuit diagrams, component part lists, descriptions, calibration instructions, and other information will assist the user's appropriately qualified technical personnel in repair. Those parts of equipment that are designated repairable can be supplied.
10. It is recommended for the performance to be checked every 2 years.
11. Please keep the cuff clean. If the cuff becomes dirty, remove it from the monitor and clean it by hand in a mild detergent, then rinse it thoroughly in cold water. Never dry the cuff in clothes dryer or iron it. Cleaning the cuff after the usage of every 200 times is recommended.
12. The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or 5 years.
13. The cuff's integrity is maintained after 1,000 open-close cycles of the closure.

PRODUCT SPECIFICATIONS

Name and Item Number	Premium Talking Digital Wrist Blood Pressure Monitor, 04-895-001
Display System	Digital display/LCD
Measuring Method	Oscillometric Method Automatic Air Inflation and Measurement
Power Source	2 × 1.5V  SIZE AAA
Measuring Range	Pressure: 0 to 300 mmHg Pulse: 40-180 beats/minute
Accuracy	Pressure: ±3 mmHg; Pulse: ±5% of reading
Diastolic Range	40-199 mmHg
Systolic Range	60-260 mmHg
Pressurization	Automatic pressurization by pump
Deflation	Automatic pressure release valve
Memory	Built-in memory enabling display of up to 120 (60 per memory bank) measurements
Automatic Power Off	Approximately 1 minute after last button operation
Battery Life	Approximately 120 times
Operation Environment	Temperature 41°F – 104°F (5°C – 40°C) Humidity ≤90% RH max
Storage Environment	Temperature -4°F – 131°F (-20°C – 55°C) Humidity ≤90% RH max
Weight	4.6 oz. (without batteries)
Dimensions	3¼" x 2¾" x 1⅛"
Wrist Circumference Range	Suitable for 5½" – 7⅝" (14 – 19.5 cm)
Accessories	Storage case, Instruction Guidebook, two AAA alkaline batteries, Quick Start Guide

Specifications are subject to change without notice.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emission		
The 04-895-001 is intended for use in the electromagnetic environment specified below. The customer or the user of the 04-895-001 should assure that it is used in such an environment		
Emissions test	Compliance	Electromagnetic environment guidance.
RF emissions CISPR 11	Group 1	The 04-895-001 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The 04-895-001 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3		

Table 2**For all ME EQUIPMENT and ME SYSTEMS**

Guidance and manufacture's declaration - electromagnetic immunity			
The 04-895-001 is intended for use in the electromagnetic environment specified below. The customer or the user of the 04-895-001 should assure that it is used in such an environment			
IMMUNITY test	IEC 60601test level	Compliance level	Electromagnetic environment guidance.
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

Table 3**For all ME EQUIPMENT and ME SYSTEMS**

Guidance and manufacture's declaration - electromagnetic immunity			
The 04-895-001 is intended for use in the electromagnetic environment specified below. The customer or the user of the 04-895-001 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601test	Compliance level	Electromagnetic environment guidance.
Radiated RF IEC 61000- 4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the 04-895-001, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Table 3 - Continued

For all ME EQUIPMENT and ME SYSTEMS

<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80MHzto 2.5 GHz</p>	<p>3 V/m</p>	<p>Recommended separation distance:</p> <p>$d=1.2 \sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d=1.3 \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level^b in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: </p>
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

Table 3 - Continued

For all ME EQUIPMENT and ME SYSTEMS

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 04-895-001 is used exceeds the applicable RF compliance level above, the 04-895-001 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the 04-895-001.

^bOver the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the 04-895-001			
The 04-895-001 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 04-895-001 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 04-895-001 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter.		
	150 kHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=1.2 \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

Table 4 - Continued

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people and animals.

**Satisfaction
Guaranteed**

Manufactured for:

HealthSmart® International

1931 Norman Drive
Waukegan, IL 60085

Register this product at LiveHealthSmart.com

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©2015 HealthSmart® International