# *Instruction Manual*

---

## **Semi-Automatic Digital Blood Pressure Monitor**

Item # 04-610-001

Please read this guidebook  
completely before operating this unit.

**DETAILED  
GUIDEBOOK**

*English • Spanish*

## Limited Three-Year Warranty

---

The warrantor guarantees that its Digital Blood Pressure Monitor will be free from manufacturing defects under normal use for three years from the original purchase date. This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the Digital Blood Pressure Monitor. This warranty is voided if the Digital Blood Pressure Monitor product is misused or abused in any manner.

If the Digital Blood Pressure Monitor fails to operate during the first three years from the time originally purchased, please call our Customer Care Help Line at 800-622-4714 and we will try to resolve the issue as soon as possible. If the warrantor determines the unit failed to operate due to a manufacturing defect, the unit will be replaced at the option of the warrantor. A return shipping and handling fee will be applied for all returns. Replacement is the sole remedy under this limited warranty. This warranty gives you specific legal rights, which vary from state to state. As a condition of this warranty, the enclosed warranty registration card must be completed and sent to us within 10 days of the purchase date or can be filled out online at [www.mabisdmi.com](http://www.mabisdmi.com).

This warranty constitutes the warrantor's only responsibility and obligation to replace materials or components. We make no other express or implied warranties, arising by operation of law or otherwise, or any warranty of merchantability or fitness for a particular use or purpose whether or not the use or purpose has been disclosed to the warrantor in specifications, drawings or otherwise, and whether or not the warrantor's products are specifically designed and/or manufactured by the warrantor for the buyer's use or purposes, except for the limited warranty stated above. The warrantor will not be responsible for any indirect, incidental, special, consequential, or punitive damages or other loss, including, but not limited to, damage to or loss of other property or equipment and personal injuries, whether to purchaser or others. The warrantor shall in no event be liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit.

1931 Norman Drive • Waukegan, IL 60085  
Toll-Free Customer Care Help Line: 1-800-622-4714  
Monday - Friday 8:00 am - 4:30 pm CST

## Index

---

Introduction	3
Product Identification & Caution	4
General Blood Pressure Information	5-7
Important Information Before Use	7-8
World Health Organization (WHO) Indicator	9
Irregular Heartbeat Detection	10
Battery Installation/Replacement	11
Setting the Date/Time	12
Applying Your Blood Pressure Cuff	13
Taking Your Blood Pressure Reading	14
Recalling Measurements in Memory	15
Description of Display Symbols	16
Troubleshooting	17
Care and Maintenance	18
Product Specifications	19
Blood Pressure Logs	25-26

## Introduction

---

*To achieve the maximum benefit from your blood pressure monitor, we recommend that you first consult with your physician or trained healthcare professional.*


Thank you for purchasing a Semi-Automatic Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

The method of measurement that your Semi-Automatic Monitor uses is called the oscillometric method. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor easy to use.

*Blood pressure readings determined with this device are equivalent to measurements obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.*

*Read all of the information in the operation guide and any other literature in the box before operating the unit.*

# Product Identification

## CAUTION:

You can stop the inflation or deflation process anytime by pressing the Release button.

# General Blood Pressure Information

## What is Blood Pressure?

Blood pressure is the pressure that is exerted by blood flowing against the walls of the blood vessels throughout your body.

Your heart, which is the center of the circulatory system, provides the force for the blood to flow or circulate. When your heart contracts or "beats", the blood is forced through the blood vessels increasing the pressure. This is the highest pressure in the cycle or what is referred to as SYSTOLIC blood pressure. In between beats, your heart relaxes and your blood pressure decreases. This is referred to as DIASTOLIC blood pressure.

This complete series of events, which occurs in a single heartbeat, is known as the CARDIAC CYCLE.

Your Semi-Automatic Blood Pressure Monitor will automatically read your blood pressure and display both systolic and diastolic readings on the screen. Your systolic will be positioned as the upper number and the diastolic reading will be the lower number.

(systolic) 120/80 (diastolic)

Blood pressure is measured in millimeters (mm) of mercury (Hg) and is generally recorded with the systolic pressure (120) listed first and the diastolic pressure (80) listed second. The numbers are typically separated by a slash mark (/) as shown.

*Both pressure readings, the SYSTOLIC and DIASTOLIC, are necessary for a physician to evaluate the status of a patient's blood pressure.*

*Please contact your physician for specific information regarding your own blood pressure.*

## What Influences Blood Pressure?

Many factors such as genetics, age, sex, altitude, physical activity, anxiety, muscular development, certain medications or even the time of day can influence blood pressure. Influences such as sleep or relaxation decrease blood pressure, while anxiety or exercise increase blood pressure.

### Why Monitor Your Blood Pressure at Home?

A visit to a physician's office can be a stressful situation for a patient. Anxiety is a known factor in raising blood pressure. This temporary occurrence of elevated blood pressure at the physician's office is commonly referred to as "white coat syndrome".

Whether or not you experience "white coat syndrome", home blood pressure monitoring provides you with the opportunity to supplement your physician's office measurements. These home readings, when taken over a period of time, can show an accurate indication of change. Furthermore, your records can assist your physician in evaluating your health and in making important decisions in the diagnosis and treatment of your condition. Because of this, it is important to take consistent, daily measurements of your blood pressure.

*The variations in your individual readings should only be interpreted by your physician or a trained health care professional.*

### WHO Blood Pressure Classifications


Standards for assessment of high or low blood pressure have been established by the World Health Organization (WHO) as shown on the following chart:

CATEGORY	SYSTOLIC (mmHg)	DIASTOLIC (mmHg)	COLOR INDICATOR
Optimal	<120	<80	GREEN
Normal	120-129	80-84	GREEN
High - Normal	130-139	85-89	GREEN
Stage 1 - Mild Borderline Hypertension	140-159	90-99	YELLOW
Stage 2 - Moderate Hypertension	160-179	100-109	ORANGE
Stage 3 - Severe Hypertension	≥180	≥110	RED
Isolated Systolic Hypertension	≥140	<90	YELLOW
	160-179	<90	ORANGE
	≥180	<90	RED

*This chart is only a general guideline. Contact your physician or trained healthcare professional to determine your NORMAL blood pressure.*

### Variations in Blood Pressure

Blood pressure is influenced by many factors and can change from moment to moment. Normally, blood pressure is lowest during sleeping periods and rises during the day. The graph below represents variations in blood pressure shown over a day with measurements taken every 5 minutes.

The dotted line represents the sleep period. The rise in blood pressure at 1 p.m. (A in the graph) corresponds to a stressful occurrence and at 7 p.m. (B in the graph) a period of exercise.

## Important Information Before Use


- Blood pressure readings should be interpreted by a physician or trained healthcare professional, who is familiar with your medical history.
- Perform your measurement in a quiet place. You should be seated and relaxed.
- Avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading. If you are exhibiting signs of stress, avoid taking your measurement until the feeling subsides.
- Rest 15 minutes prior to taking a reading.
- Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
- Keep the monitor stable during measurements. Remain still; do not talk during measurements.
- Record your daily readings on the chart provided.
- Take your readings at the same time, each day (or as recommended by your physician).
- Wait a minimum of 15 minutes between readings. The wait time may vary depending on individual physiological characteristics.
- This device is intended for adult use.
- The inflation or deflation process can be stopped by pressing the Release button.
- People with serious arrhythmia are discouraged from using this Electronic Sphygmomanometer.


## Important Information Before Use

13. Use only the cuff supplied by the manufacturer to prevent measurement errors caused by biochemical hazard from occurring.
14. To avoid infection, do not share cuff with other users.
15. This Electronic Sphygmomanometer should never be used on infants or young children. Consult your physician before use on older children.
16. Prolonged over-inflation (cuff pressure exceeds 300 mmHg or maintained above 15 mmHg for longer than 3 minutes) of bladder may cause ecchymoma of the arm.
17. Avoid application of cuff over a wound or inflammatory diseases.
18. Avoid application of cuff on any limb where intravascular access or therapy, or an arteriovenous (A-V) shunt, is present.
19. Avoid application of the cuff on the arm on the side of a mastectomy.
20. May not be simultaneously used with other monitoring medical equipment on the same limb;
21. Need to check the blood circulation of the user.
22. The monitor might not meet performance specifications or may cause safety hazard if stored or used outside the specified temperature and humidity ranges.
23. For information regarding potential electromagnetic or other interference between the blood pressure monitor and other devices together with advice regarding avoidance of such interference, please see part ELECTROMAGNETIC COMPATIBILITY INFORMATION.
24. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
  - Reorient or relocate the receiving antenna.
  - Increase the separation between the equipment and receiver.
  - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
  - Consult the dealer or an experienced radio/TV technician for help.

## World Health Organization (WHO) Indicator

The WHO Indicator, located on the left side of the display, will identify your reading within a range in accordance with the guidelines set by the World Health Organization (WHO). See the chart below for reference.

SYS  DIA 	<b>Stage 3 Hypertension</b> Systolic pressure: $\geq 180$ mmHg -or- Diastolic pressure: $\geq 110$ mmHg
SYS  DIA 	<b>Stage 2 Hypertension</b> Systolic pressure: 160-179 mmHg -or- Diastolic pressure: 100-109 mmHg
SYS  DIA 	<b>Stage 1 Hypertension</b> Systolic pressure: 140-159 mmHg -or- Diastolic pressure: 90-99 mmHg
SYS  DIA 	<b>High Normal</b> Systolic pressure: 130-139 mmHg -or- Diastolic pressure: 85-89 mmHg
SYS  DIA 	<b>Normal</b> Systolic pressure: 120-129 mmHg -or- Diastolic pressure: 80-84 mmHg
SYS  DIA 	<b>Optimal</b> Systolic pressure: $< 120$ mmHg -and- Diastolic pressure: $< 80$ mmHg

## Irregular Heartbeat Detection

Your digital blood pressure monitor features Irregular Heartbeat Detection. This feature allows users to accurately monitor blood pressure even if an irregular heartbeat should occur.

When an irregular heartbeat is detected, the IHB icon will appear on the display, Fig. 1.

**NOTE:** Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat if this symbol appears frequently at the end of a reading.

If Irregular Heartbeat (IHB) brought by common arrhythmias is detected during blood pressure measurement, a signal will be displayed. Under this condition, the device may continue to function, but the results may not be accurate. It's suggested that you consult with your physician for accurate assessment.


Fig. 1

## Battery Installation/Replacement

1. Press down on tab and lift the cover.
2. Insert or replace 4 "AA" batteries into the compartment, matching the indicated polarity symbols, Fig. 1.
3. To close, match up and connect the bottom hooks of the cover, then push the top end.

The monitor, the batteries and the cuff must be disposed of according to local regulations at the end of their usage.


Fig. 1

**NOTE:** Replace the batteries when the Low Battery symbol appears on the display or when the display remains blank when the POWER button is pressed.

4. It is recommended to remove the batteries if the unit will not be used for an extended period of time.

## Setting the Date/Time

**NOTE:** When the unit is off, the date and time will still appear at the top of the screen.

1. When the unit is off, press and release both the POWER and MEMORY button at the same time. The month will begin to flash, Fig. 1.
2. Press the MEMORY button to change until the desired number appears.

**NOTE:** To scroll through the numbers quicker, hold down the MEMORY button.

3. Press the POWER button to set the entry.
4. The date, Fig. 2, followed by the hour and minutes will blink. Repeat steps 2 and 3 to adjust the date/time setting.
5. After setting the minutes, the unit will automatically exit out of the date/time setting mode, Fig. 3.

**NOTE:** The time is expressed in 24hr military time. The chart below explains the time conversion.

Regular Time (12 hours)	Military Time (24 hours)	Regular Time (12 hours)	Military Time (24 hours)
12 am Midnight	0	12 pm Noon	12
1 am	1	1 pm	13
2 am	2	2 pm	14
3 am	3	3 pm	15
4 am	4	4 pm	16
5 am	5	5 pm	17
6 am	6	6 pm	18
7 am	7	7 pm	19
8 am	8	8 pm	20
9 am	9	9 pm	21
10 am	10	10 pm	22
11 am	11	11 pm	23


Fig. 1

Fig. 2


Fig. 3

## Applying Your Blood Pressure Cuff

Before applying your blood pressure cuff, be sure you have selected the appropriate cuff:

**Standard Cuff** fits arm circumference 11-3/4" – 16-1/2" (30-42 cm)

**Large Cuff** fits arm circumference 16-1/2" – 18-7/8" (42-48 cm)

**NOTE:** If for any reason you are unable to or should not use your left arm, please modify the following instructions and apply the cuff to your right arm. Your physician can tell you which arm is best for you to use.

1. Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
2. Be seated with your feet flat on the floor.
3. Position the cuff on a solid surface with the tubing facing UP and AWAY from you, Fig. 1. The metal bar on the cuff should be to the left of the tubing.
4. Widen or open the cuff by pulling the bottom of the cuff towards the right, Fig. 2. This should open the cuff, creating a cylinder. Do not extend the cuff beyond the metal bar.
5. Insert your arm into the cuff (cylinder). Position the (Φ) mark over the main artery (on the inside of your arm), Fig. 3.
6. The bottom edge of the cuff should be positioned approximately one inch ABOVE the elbow joint, Fig. 4.
7. Reaching underneath your left arm with your right hand, pull the end of the cuff towards your body to tighten the cuff, Fig. 5. Wrap and secure the cuff making sure that the (Φ) mark remains as shown, Fig. 3.

Fig. 1

Fig. 2


Fig. 3


Fig. 4


Fig. 5

**NOTE:** The cuff size is suitable for use when the vertical "Index" mark is within the horizontal "OK" range. A different sized cuff is needed if the "Index" mark is outside of the "OK" range.

Proceed only after reading the previous sections of this manual.

## Taking Your Blood Pressure Reading

1. Position the monitor on a flat, stable surface with the digital display panel in view. Insert the cuff tubing connector into the left side of the monitor, *Fig. 1*. Avoid compression or restriction of the connection tubing during measurement which may cause inflation error or harmful injury due to continuous cuff pressure.


Fig. 1

2. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, *Fig. 2*. Relax your left hand, sit comfortably and don't cross your legs. The middle of the cuff should be at the level of the right atrium of the heart.


Fig. 2

3. Press the POWER Button again to turn the unit on.
4. This unit will run a self-test, *Fig. 3*, and show the last measured reading, *Fig. 4*, or if there is no previous measurement, "0" as show in *Fig. 5*.

Fig. 3


Fig. 4


Fig. 5

5. When the display appears as shown, *Fig. 6*, the monitor is ready.

Note: If a downward arrow (↓) flashes, there is air remaining in the cuff. The air will be released by the "Automatic Electronic Air Release Valve" until "0" appears on the display.

To release the air quickly press the release button.

6. Using your right hand, squeeze and release the inflation bulb at a steady pace. The inflation pressure must be 50-60 mmHg above your normal systolic pressure. If the pressure in the cuff is insufficient, an upward arrow (↑) will flash. Continue to steadily inflate the cuff to the proper level. Once the unit beeps, stop inflating the cuff. The unit will automatically begin to deflate.


Fig. 6

7. Upon completion, your blood pressure (systolic and diastolic) measurement will display with the downward arrow (↓) flashing to let you know to press the release button to release the residual air quickly, *Fig. 7*.

Fig. 7

8. Your measurement will be automatically saved in the memory.
9. To conserve the batteries, press the POWER button to turn the power off. Otherwise, the unit will automatically shut off after approximately 1 minute.
10. Disconnect the tubing from the monitor prior to storing.

## Recalling Measurements in Memory

You can recall up to 60 total measurements in memory. These measurements can be shared with your physician or trained healthcare professional.

1. Press and release the MEMORY button. It will show how many measurements are stored within the bank, followed by the last result, *Fig. 1*.
2. Press the MEMORY button again to scroll through the saved measurements. After you've scrolled through all the readings, the unit will automatically turn off.


Fig. 1

**NOTE:** When the number of readings exceeds 60, the oldest data will be replaced with the new record.

**To erase all stored measurements:**

1. Press and release the MEMORY button. The unit will display how many measurements are stored within the bank.
2. Press and hold down the MEMORY button until two horizontal lines flash. This indicates that all measurements have been erased, *Fig. 2*.


Fig. 2

## Description of Display Symbols

SYS 	Indicates a measurement within Stage 3 Hypertension.	 Shows the pulse rate per minute.
DIA 		 Occurs when an error was made during measurement.
SYS 	Indicates a measurement within Stage 2 Hypertension.	 Appears when there is an irregular heartbeat detected during measurement.
DIA 		 Appears when the cuff is unstable or it has too much air left in it. Once the arrow stops flashing, the unit is ready to take a measurement.
SYS 	Indicates a measurement within Stage 1 Hypertension.	 Appears when the unit is ready or a measurement error has occurred. Restart the measurement if an error occurs.
DIA 		
 Appears when batteries should be replaced.		 Appears when recalling a measurement from memory.

## Explanation of Symbols

Symbol for "THE OPERATION GUIDE MUST BE READ" (The sign background color: blue. The sign graphical symbol: white)

Symbol for "WARNING"

Symbol for "TYPE BF APPLIED PARTS" (The cuff is type BF applied part)

Symbol for "ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local Authority or retailer for recycling advice".

Symbol for "KEEP DRY"

## Troubleshooting

If any abnormality occurs during use, please check and correct the following:

Condition	Correction
Display is blank when power is on.	Check and correct the polarity of the installed batteries. Reinstall or replace batteries.
Measurement incomplete or abnormally low or high values displayed.	Review and follow "Applying Your Blood Pressure Cuff" and "Taking Your Blood Pressure Reading" sections. Replace the batteries.
EE symbol displays.	Do not talk or move during the measurement. Review "Applying Your Blood Pressure Cuff" section.
Measurements are different from those typically measured by physician or every measurement reading is different.	Remember, blood pressure readings are influenced by physical and mental conditions and/or even the time of day. Daily measurements should be interpreted by your physician.

## Care and Maintenance

---

## Product Specifications

---

<b>Name and model number</b>	Semi-Automatic Standard Digital Blood Pressure Arm Monitor, 04-610-001
<b>Display system</b>	Digital display/LCD
<b>Measuring method</b>	Oscillometric
<b>Power source</b>	4 x 1.5V  Size AA
<b>Measuring range</b>	Pressure: 0-300 mmHg Pulse: 40-180 beats/minute
<b>Accuracy</b>	Pressure: $\pm 3$ mmHg Pulse: $\pm 5\%$ of reading
<b>Diastolic</b>	40-199mmHg
<b>Systolic</b>	60-260mmHg
<b>Inflation</b>	Manual Inflation
<b>Deflation</b>	Automatic Pressure Release Valve/Manual fast release
<b>Memory</b>	Built-in memory enabling display of up to 30 measurements
<b>Automatic shut-off</b>	Approximately 1 minute after last button operation
<b>Battery Life</b>	Approximately 180 uses
<b>Operation environment</b>	Temperature 41°F – 104°F (5°C – 40°C) Humidity $\leq 90\%$
<b>Storage environment</b>	Temperature -4°F – 131°F (-20°C – 55°C) Humidity $\leq 90\%$
<b>Monitor dimensions</b>	5" x 5-1/8" x 2-3/8"
<b>Weight</b>	11.11 oz. (without batteries)
<b>Arm circumference range</b>	Adult size cuff fits arm circumference: 11-3/4" – 16-1/2" (30-42 cm)
<b>Accessories</b>	Cuff, detailed guidebook, quick start guide, blood pressure log

*Specifications are subject to change without notice*

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 1  
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions		
The 04-610-001 is intended for use in the electromagnetic environment specified below. The customer or the user of the 04-610-001 should assure that it is used in such an environment		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The 04-610-001 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The 04-610-001 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions	Not Applicable	
IEC 61000-3-2		
Voltage fluctuations/flicker emissions IEC 61000-3-3		

Table 2  
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic immunity			
The 04-610-001 is intended for use in the electromagnetic environment specified below. The customer or the user of the 04-610-001 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U <sub>T</sub> is the a.c. mains voltage prior to application of the test level.			

**Table 3**  
For all ME EQUIPMENT and ME SYSTEMS

<b>Guidance and manufacture's declaration - electromagnetic immunity</b>			
The 04-610-001 is intended for use in the electromagnetic environment specified below. The customer or the user of the 04-610-001 should assure that it is used in such an environment.			
<b>IMMUNITY test</b>	<b>IEC 60601 test</b>	<b>Compliance level</b>	<b>Electromagnetic environment-guidance</b>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the 04-610-001, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

**Table 3 - Continued**  
For all ME EQUIPMENT and ME SYSTEMS

Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<b>Recommended separation distance:</b>
			$d=1.2 \sqrt{P}$ 80 MHz to 800 MHz  $d=1.3 \sqrt{P}$ 800 MHz to 2.5 GHz  Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level <sup>b</sup> in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the 04-610-001.			
The 04-610-001 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 04-610-001 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 04-875-001 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter.		
	150 kHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=1.2 \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people and animals.			

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1931 Norman Drive South • Waukegan, IL 60085  
Toll-Free Customer Care Help Line: 1-800-622-4714  
Monday - Friday 8:00 am - 4:30 pm CST

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